
EVAPORADOR DE PELICULA ASCENDENTE

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QTS 4 Evaporation - Climbing Film and Natural Circulation Evaporator

QVF Process Systems, Inc. 224 North Main Street, Building F, Section 4, Horseheads, New York 14845. Telephone: (607) 796-4800, Fax: (607) 739-6717, E-Mail:info@qvfps.com. <http://www.qvfps.com/qts4.html>

Specification

A QVF evaporation teaching system module comprising a 3000mm x 25mm diameter long tube evaporator, cyclone separator, vapor condenser, product and condensate receivers, all manufactured in borosilicate glass. Designed to operate as either a climbing film or natural circulation evaporator for atmospheric or vacuum operation. Supplied with all valves, controllers, instrumentation, service lines and supporting structure. Complete with an installation and maintenance manual and tool kit.

Experimental Capabilities

- Mass and heat balances
- Determination of evaporative heat transfer coefficients
- Influence of vacuum levels on evaporation rates
- Influence of flowrate on forced circulation heat transfer coefficients
- Comparison of climbing film, natural and forced circulation evaporation

Water is used in the experiments. Climbing film evaporation is particularly suited to the concentration of heat sensitive materials such as milk, fruit juices, blood plasma, etc. and for desalination studies.

Description

The plant is fitted with a steam heated long tube calandria evaporator. Feed material from the feed vessel enters at the calandria base, its flowrate being determined by means of an in-line variable area flow meter, and is partially evaporated. Vapor and liquid pass into the cyclone which separates the phases. Vapor passes to a water cooled condenser, the condensate being collected in a receiver which allows condensate removal while operating under vacuum.

When operating as a climbing film evaporator, feed material partially evaporates as the liquid film is seen to climb the tube. Concentrate from the cyclone is collected in a graduated receiver as final product or is recirculated back to the calandria

In natural circulation mode, the evaporator is run flooded and concentrate continually passes back to the evaporator, evaporated product being replaced intermittently from the feed vessel.

The optional equipment consisting of a pump, variable area flow meter and interconnecting pipework allows forced circulation evaporation to be studied over a wide range of flowrates.

Plant operation is controlled from a free-standing control panel fitted with a mimic diagram, flow controllers, temperature and pressure indicators and services controllers.

All major process items and pipework are manufactured from QVF borosilicate glass with other contact parts in PVDF, PTFE and stainless steel.

Options: Either steam or electric heating.

Electrical System: The electrical system can be designed to meet your local requirements.

Equipment Specification: Calandria: 3000mm x 25mm diameter
Condenser: 150mm dia x 1.5m²
Concentrate receive (1): 7 L graduated in 50 ml increments
Product receivers (2) 5L

Instrumentation: Water: Control valve, pressure gauge
Vacuum: Control valve, pressure gauge
Steam: Reducer, control valve, pressure gauge
Process: Two temperature indicators

Utilities Requirements: Steam: 45 lb/hr at 2.1 barg.
Water: 1250 L/h at 1 barg
Vacuum: 100 torr

Space Requirements : (height x floor area) 5.2 x 2.0 1.0 meters.